

FLAME PHOTOMETER FP8000 SERIES

DETERMINING THE CONCENTRATION OF ALKALINE AND EARTH ALKALINE ELEMENTS IN THE LABORATORY AND PROCESS ANALYSIS



FLAME PHOTOMETER FP8000 SERIES



The A.KRÜSS flame photometers allow simple and cost-effective determination of alkali and alkaline earth element concentrations in aqueous solutions. They are mainly used in laboratory and process operations.

With the FP8000 series A.KRÜSS Optronic presents worldwide unique measuring devices for simultaneous measurement of the elements Na, K, Ca and Li. Five measurement channels are available to be individually configured. Thanks to the proven measurement procedure as well as to modern device electronics and innovative process control, our flame photometers achieve highly accurate and reproducible results. The devices can cope with about 300 measurements per hour and thus clearly demonstrate the international top position!

KEY FEATURES

- Economic, robust and durable solution for use in laboratory and process operations
- Maximum operational safety through intelligent security mechanisms
- Flame photometry simultaneous measurement of up to five alkaline and alkaline earth elements
- Fast measurement: up to 300 measurements per hour international top position!
- User administration with two authorisation levels
- Traceability of all measured data and device-specific data
- Complete traceability of measurement results
- Extensive interfacing options and comfortable handover of measurement results
- Compliance with international norms and standards such as GMP / GLP and 21 CFR Part 11

SOLUTIONS FOR VARIOUS NEEDS



	FP8400	FP8500	FP8600	FP8700
APPLICATIONS	Standard version for Laboratory Economic entry solution	Version for Process Remote control of multiple devices via centre process computer	Version with automated sample feed	Version with automated sample feed and dilution
	SRUSS SR	XXXX		
SAMPLER			✓	√
DILUTER				✓
MIXER				√
VALVE BOX		✓	✓	

EQUIPMENT AND FIELDS OF APPLICATION



	FP8400	FP8500	FP8600	FP8700	
DETERMINING CONCENTRATION	Determining the concentrations of alkaline and alkaline earth elements in aqueous solutions for laboratory and process operations Na, K, Ca, Li - Propane or acetylene				
EQUIPMENT	Basic device	Basic device+	Basic device +	Basic device +	
		External valve box with control unit (Online Measurement)	Sampler and external sample table and control unit	External automation system consisting of: sampler, diluter and mixer	
FIELDS OF APPLICATION	Laboratory analysis Quality control and evidence Food laboratories Examination offices Chemical development	Process analysis Process monitoring	Quality control and evidence Food laboratories Examination offices	Quality control and evidence Food laboratories Examination offices	
	For operation with low sample throughput Very high precision	Recommended for online measurements Designed for automated continuous measurements	Automated procedures for operation with high sample throughput Samples are measured undiluted Little support from users needed High precision	Automated procedures for operation with high sample throughput Samples are measured diluted Little support from users needed High precision Small amount of sample	

EASY HANDLING AND OPERATING COMFORT



OPERATION

Our flame photometers are easy and intuitive to use on a large touchscreen display and with USB mouse.

The FP8000 series is designed for using propane or acetylene as fuel gas. Purified and dried compressed air is required for operation.

CALIBRATION, SAMPLES, CONTROLS AND STANDARDS


Aqueous samples, controls and standards can be measured manually or automatically. Calibration is carried out via the standards.

Measurement results can be checked with the controls.

For reviewing, emission and concentration values of calibration are presented in graphed and tabulated order.

USER ADMINISTRATION

All devices feature a user administration with two authorization levels, which can selectively be activated or deactivated.



MEASUREMENT PRINCIPLE, FUNCTIONALITY, INNOVATIVE PROCESS CONTROL



MEASUREMENT PRINCIPLE OF FLAME PHOTOMETERS

The flame photometers of the FP8000 series apply a relative measurement method according to the principle of flame atomic emission spectroscopy (F-AES).

FUNCTIONALITY OF THE FP8000 SERIES

The FP-8000 devices come with five independent measurement channels and another, separate channel for monitoring the flame. The five measurement channels can be fitted with appropriate filters. The devices do not necessarily need a guideline, measurement can be done both with and without reference element.

INNOVATIVE PROCESS CONTROL

The flame photometers can be used as table top device or be integrated into a process environment via analogue current interface.

Several locally separated units can be controlled via one centre process computer.

The units are designed for 24-hour operation.



FP8500

HIGH SAFETY STANDARD FOR MAXIMUM OPERATIONAL SAFETY



Various safety mechanisms ensure particularly high operational safety. Such as a safety circuit for interrupting the supply of gas preventing leakage of unburned gases, or as a flame recoil blocker preventing the backward burning of an extinguished flame.

SOME ESSENTIAL SAFETY MECHANISMS

Automated ignition

• Controlled and safe ignition sequence

Optical infrared flame monitoring

• Continuous monitoring of flame via infrared sensor

Continuous monitoring of gas pressure

• Fuel, Oxidant and Cooling air



HIGH SAFETY STANDARD FOR MAXIMUM OPERATIONAL SAFETY



SOME ESSENTIAL SAFETY MECHANISMS

Automated gas shut down

- In case of an extinguished flame
- In case of regular system shut down
- In case of power breakdown
- In case of low fuel-gas, oxidant or cooling air pressure
- The monitoring of the safety functions is performed by an independent safety circuit

Exhaust emission

 \bullet To avoid injury due to hot exhaust gases, exhaust gas temperatures are held at a cooling air temperature below ${\sim}50\,^{\circ}{\rm C}$



INTERFACES AND DATA EXPORT



RS-232	USB	ethernet
Serial ASCII printer	USB flash drive Barcode scanner Keyboard Mouse	Network printer Network share LIMS

- Print-out on serial ASCII printer at RS-232 interface
- Print-out on network printer in PDF or GS format
- Print-out as PDF on USB flash drive or to network share
- Export of the measurement results in HTML or CSV format on
- USB flash drive or to network share
- Easy integration into existing networks (DHCP client) or a LIMS

SERVICE AND MAINTENANCE - YOUR BENEFITS AS AN A.KRÜSS CUSTOMER



- On site services by A.KRÜSS Optronic or certified service partner: :
 - o Service
 - Maintenance
 - O Initial installation and commissioning
 - o DQ/IQ/OQ/PQ
 - O Calibration and adjustment
- Calibration and adjustment with defined standard solutions
- Maintenance contracts, calibration certificates
- Individual training and application consulting on site
- Customer-specific customisation of devices
 (e. g. optimisation of measuring range)
- Spare parts and accessories directly from the manufacturer



FEATURES AND SPECIFICATIONS



- Economic, robust and durable solution for use in laboratory and process operation
- Flame photometry simultaneous measurement of up to five alkaline and alkaline earth elements
- Use as table top unit or integration into a process environment (analogue current interface)
- TFT-Display with integrated touchscreen, 8,4" (21 cm), 800x600 Pixel
- Easy and intuitive operation, available in German and English language (others on request)
- Aqueous samples, controls and standards are measured manually or automatically
- Graphic and tabular presentation of emission and calibration curves
- Method Management: create an arbitrary number of methods to enable the respectively desired parameters to each sample
- Available method parameters: elements, concentration, calibration, controls, comment
- Fast measurement: up to 300 measurements per hour
- User administration: optional setup, two authorisation levels can be activated or deactivated depending on the particular requirements
- Integrated measured data storage: document all measured values of the last 999 measurements in circular buffer
- Seamless documentation: complete recording of all measured data as well as system or method settings
- Spare parts and accessories directly from the manufacturer
- Measurement principle: F-AES flame atomic emission spectroscopy
- Maximum operational safety through intelligent security mechanisms

TECHNICAL DATA IN ppm



		FP8400	FP8500	FP8600	FP8700
measuring range	Na	0.01-4500 ppm	0.01-4500 ppm	0.01-4500 ppm	0.1-45000 ppm
	K	0.01-4500 ppm	0.01-4500 ppm	0.01-4500 ppm	0.1-45000 ppm
	Li	0.01-4500 ppm	0.01-4500 ppm	0.01-45000 ppm	0.1-45000 ppm
	Ca	0.50-4500 ppm	0.50-4500 ppm	0.50-45000 ppm	5.0-45000 ppm
DETECTION	Na	0.01 ppm	0.01 ppm	0.01 ppm	0.1 ppm
LIMIT	K	0.01 ppm	0.01 ppm	0.01 ppm	0.1 ppm
	Li	0.01 ppm	0.01 ppm	0.01 ppm	0.1 ppm
	Ca	0.03 ppm	0.03 ppm	0.03 ppm	0.3 ppm
PRECISION	Na	0.2 % at 40 ppm	0.2 % at 40 ppm	0.4 % at 40 ppm	0.6 % at 40 ppm
	K	0.2 % at 40 ppm	0.2 % at 40 ppm	0.4 % at 40 ppm	0.6 % at 40 ppm
	Li	0.2 % at 40 ppm	0.2 % at 40 ppm	0.4 % at 40 ppm	0.6 % at 40 ppm
	Ca	0.2 % at 40 ppm	0.2 % at 40 ppm	0.4 % at 40 ppm	0.6 % at 40 ppm
ACCURACY	Na	1 % at 40 ppm			
	K	1 % at 40 ppm			
	Li	1 % at 40 ppm			
	Ca	1 % at 40 ppm			

TECHNICAL DATA IN mmol/l



		FP8400	FP8500	FP8600	FP8700
measuring range	Na	0.0004-200 mmol/l	0.0004-200 mmol/l	0.0004-200 mmol/l	0.004-2000 mmol/l
	K	0.0003-110 mmol/l	0.0003-110 mmol/l	0.0003-110 mmol/l	0.003-1100 mmol/l
	Li	0.0014-600 mmol/l	0.0014-600 mmol/l	0.0014-600 mmol/l	0.014-6000 mmol/l
	Ca	0.0125-110 mmol/l	0.0125-110 mmol/l	0.0125-110 mmol/l	0.125-110 mmol/l
DETECTION	Na	0.0004 mmol/l	0.0004 mmol/l	0.0004 mmol/l	0.0004 mmol/l
LIMIT	K	0.0003 mmol/l	0.0003 mmol/l	0.0003 mmol/l	0.0003 mmol/l
	Li	0.0014 mmol/l	0.0014 mmol/l	0.0014 mmol/l	0.0014 mmol/l
	Ca	0.0075 mmol/l	0.0075 mmol/l	0.0075 mmol/l	0.0075 mmol/l
PRECISION	Na	0.2% at 1.74 mmol/l	0.2 % at 1.74 mmol/l	0.4 % at 1.74 mmol/l	0.6 % at 1.74 mmol/l
	K	0.2% at 1.03 mmol/l	0.2 % at 1.03 mmol/l	0.4 % at 1.03 mmol/l	0.6 % at 1.03 mmol/l
	Li	0.2% at 5.71 mmol/l	0.2 % at 5.71 mmol/l	0.4 % at 5.71 mmol/l	0.6 % at 5.71 mmol/l
	Ca	0.2% at 1.00 mmol/l	0.2 % at 1.00 mmol/l	0.4 % at 1.00 mmol/l	0.6 % at 1.00 mmol/l
ACCURACY	Na	1% at 1.74 mmol/l	1 % at 1.74 mmol/l	1 % at 1.74 mmol/l	1 % at 1.74 mmol/l
	K	1% at 1.03 mmol/l	1 % at 1.03 mmol/l	1 % at 1.03 mmol/l	1 % at 1.03 mmol/l
	Li	1% at 5.71 mmol/l	1 % at 5.71 mmol/l	1 % at 5.71 mmol/l	1 % at 5.71 mmol/l
	Ca	1% at 1.00 mmol/l	1 % at 1.00 mmol/l	1 % at 100 mmol/l	1 % at 1.00 mmol/l

TECHNICAL DATA



	FP8400	FP8500	FP8600	FP8700	
CALIBRATION	Linear with 2 standards				
		Non linear with 6-8 standa	rds, cubic approximation		
DRIFT					
REFERENCE		Lithium guide line	e with 5 mmol/l		
SAMPLE VOLUME	2,5 ml	2,5 ml	2,5 ml	0,25 ml	
COMBUSTION GAS	Propane: recommended for alkaline elements				
		Butane: possible sub	stitute for propane		
	Acetylene for flame photometry: recommended for alkaline earth elements			ements	
DISPLAY	TFT-	Display with integrated 8,4" T	TFT touchscreen, 800x600 Pix	cel	
INTERFACES	2 x USB flash drive 1 x Ethernet 1 x RS-232 for printer	2 x USB flash drive 1 x Ethernet 1 x RS-232 for printer Upgrade 4-20 mA analogue, passive current interface	2 x USB flash drive 1 x Ethernet 1 x RS-232 for printer	2 x USB flash drive 1 x Ethernet 1 x RS-232 for printer	
ACCESSORIES	24 character plain paper printer				
	Pressure reducing station with filter unit				
	Various standards and system solutions				
	Spare parts and accessories directly from the manufacturer				
UPGRADEABILITY	Retrofit kits to other models available				

COMPLIANCE WITH INTERNATIONAL NORMS AND STANDARDS



21 CFR Part 11

GMP / GLP

	SOME SELECTED DIN STANDARDS
DIN EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use Part 1: General requirements
DIN EN 61010-2-061	Safety requirements for electrical equipment for measurement, control, and laboratory use. Particular requirements for laboratory atomic spectrometers with thermal atomization + ionization
DIN EN 61000-3-2	Electromagnetic compatibility (EMC) — Limits for harmonic current emissions (Equipment input current)
DIN EN 61000-4-2	Electromagnetic compatibility (EMC) Testing and measurement techniques Electrostatic discharge immunity test
DIN EN 61000-4-4	Electromagnetic compatibility (EMC) Testing and measurement techniques Electrical fast transient / burst immunity test
DIN EN 61000-4-5	Electromagnetic compatibility (EMC) Testing and measurement techniques Surge immunity test
DIN EN 50082-1	Electromagnetic compatibility (EMC) Generic immunity standard Residential, commercial and light industry
DIN EN 55011	Industrial, scientific and medical equipment Radio-frequency disturbance characteristics – Limits and methods of measurement



Picture credits: S. 7 und 8: Flame, © jihane37 #80339975/www.fotolia.de S. 10: © gstockstudio/www.fotolia.de

A.KRÜSS Optronic GmbH Alsterdorfer Straße 276-278 22297 Hamburg | Germany

Tel +49 40 514317-0

Fax +49 40 514317-60

E-Mail info@kruess.com

Web www.kruess.com

