

PT 31 Electronic Peltier Thermostat



PT 31

PT 31 Adjustable Electronic Peltier Thermostat

This electronic water-bath thermostat with Peltier element is a versatile, high-performance instrument. In one application, for example, it can be used to set the correct refractometer temperature. It is extremely robust, compact and easy to operate. Thanks to its small footprint it does not take up valuable space in the laboratory.

- ▶ Temperature range : 8° - 40 °C, adjustable
- ▶ Accuracy of temperature: ± 0.2°C
- ▶ Resolution: 0.1°C
- ▶ Pump performance: 20 l/h
- ▶ Max.pump pressure: 2000 Pa
- ▶ Cooling power: 15 W
- ▶ Heating power: 30 W
- ▶ Power supply: 115-230 V AC
- ▶ Dimension L/W/H: 140x80x210 mm
- ▶ Weight: 1.5 kg
- ▶ With pump and hose connections



PT 31



A. Krüss Optronic GmbH

Alsterdorferstraße 220 · 22297 Hamburg · Germany
 Telefon +49 40 514 31 70 · Telefax +49 40 51 25 22
 Internet: www.kruess.com · E-Mail: info@kruess.com